


Product Datasheet

Mouse Cardiotrophin-1 ELISA Kit (orb1473388)

| | |
|----------------------------|---|
| Description | Mouse Cardiotrophin-1 ELISA Kit |
| Reactivity | Mouse |
| Range | 156 pg/ml - 10000 pg/ml |
| Tested Applications | ELISA |
| Target | CTF1 |
| Storage | Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles. |
| Note | For research use only |
| Application notes | Detection Method: Colorimetric, Sample Types: Cell Culture Supernates, Serum, Plasma |
| Sample Types | Cell Culture Supernates, Serum, Plasma |
| Uniprot ID | Q60753 |
| Sensitivity | 80 pg/ml |
| Entrez | 13019 |
| Expiration Date | Please enquire. |



Representative Standard Curve