


Product Datasheet

Mouse ARMETL1 ELISA Kit (orb1473396)

Description	Mouse ARMETL1 ELISA Kit
Reactivity	Mouse
Range	15.6 pg/ml - 1000 pg/ml
Tested Applications	ELISA
Target	CDNF
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Detection Method: Colorimetric, Sample Types: Cell Culture Supernates, Serum, Plasma
Sample Types	Cell Culture Supernates, Serum, Plasma
Uniprot ID	Q8CC36
Sensitivity	8 pg/ml
Entrez	227526
Expiration Date	Please enquire.



Representative Standard Curve