


## Product Datasheet

### Human PLA2G2A ELISA Kit (orb1473433)

|                            |   |
|----------------------------|---|
| <b>Description</b>         | Human PLA2G2A ELISA Kit   |
| <b>Reactivity</b>          | Human   |
| <b>Range</b>               | 156 pg/ml - 10000 pg/ml   |
| <b>Tested Applications</b> | ELISA   |
| <b>Target</b>              | PLA2G2A   |
| <b>Storage</b>             | Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles. |
| <b>Note</b>                | For research use only   |
| <b>Application notes</b>   | Detection Method: Colorimetric, Sample Types: Cell Culture Supernates, Serum, Plasma          |
| <b>Sample Types</b>        | Cell Culture Supernates, Serum, Plasma  |
| <b>Uniprot ID</b>          | <b>P14555</b>   |
| <b>Sensitivity</b>         | 80 pg/ml  |
| <b>Entrez</b>              | <b>5320</b>   |
| <b>Expiration Date</b>     | Please enquire.   |



Representative Standard Curve