


Product Datasheet

Human PDGFD ELISA Kit (orb1473435)

Description	Human PDGFD ELISA Kit
Reactivity	Human
Range	31.2 pg/ml - 2000 pg/ml
Tested Applications	ELISA
Target	PDGFD
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Detection Method: Colorimetric, Sample Types: Cell Culture Supernates, Serum, Plasma
Sample Types	Cell Culture Supernates, Serum, Plasma
Uniprot ID	Q9GZP0
Sensitivity	15 pg/ml
Entrez	80310
Expiration Date	Please enquire.



Representative Standard Curve