


Product Datasheet

Human Neutrophil Elastase ELISA Kit (orb1473440)

Description	Human Neutrophil Elastase ELISA Kit
Reactivity	Human
Range	312 pg/ml - 20000 pg/ml
Tested Applications	ELISA
Target	ELANE
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Detection Method: Colorimetric, Sample Types: Cell Culture Supernates, Serum, Plasma
Sample Types	Cell Culture Supernates, Serum, Plasma
Uniprot ID	P08246
Sensitivity	150 pg/ml
Entrez	1991
Expiration Date	Please enquire.



Representative
Standard
Curve