


Product Datasheet

Human Kallikrein 8 ELISA Kit (orb1473447)

Description	Human Kallikrein 8 ELISA Kit
Reactivity	Human
Range	156 pg/ml - 10000 pg/ml
Tested Applications	ELISA
Target	KLK8
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Detection Method: Colorimetric, Sample Types: Cell Culture Supernates, Serum, Plasma, Cell Lysates
Sample Types	Cell Culture Supernates, Serum, Plasma, Cell Lysates
Uniprot ID	O60259
Sensitivity	80 pg/ml
Entrez	11202
Expiration Date	Please enquire.



Representative Standard Curve