


Product Datasheet

Human Hepcidin ELISA Kit (orb1473462)

Description	Human Hepcidin ELISA Kit
Reactivity	Human
Range	7.8 pg/ml - 500 pg/ml
Tested Applications	ELISA
Target	HAMP
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Detection Method: Colorimetric, Sample Types: Cell Culture Supernates, Serum, Plasma, Urine
Sample Types	Cell Culture Supernates, Serum, Plasma, Urine
Uniprot ID	P81172
Sensitivity	4 pg/ml
Entrez	57817
Expiration Date	Please enquire.



Representative
Standard
Curve