


Product Datasheet

Mouse Insulin ELISA Kit (orb1473551)

Description	Mouse Insulin ELISA Kit
Reactivity	Mouse
Range	156 pg/ml - 10000 pg/ml
Tested Applications	ELISA
Target	INS
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Detection Method: Colorimetric, Sample Types: Cell Culture Supernates, Serum, Plasma, Tissue Homogenates
Sample Types	Cell Culture Supernates, Serum, Plasma, Tissue Homogenates
Uniprot ID	P01325 , P01326
Sensitivity	80 pg/ml
Entrez	16333 , 16334
Expiration Date	Please enquire.



Representative Standard Curve