


Product Datasheet

Total Collagen Microplate Assay Kit (orb1473553)

Description	Total Collagen Microplate Assay Kit
Reactivity	All
Range	0.01 mg/ml - 1 mg/ml
Tested Applications	FA
Storage	Shipped and store at 4°C for 6 months.
Note	For research use only
Application notes	Detection Method: Colorimetric, Sample Types: Urine, Serum, Plasma, Tissue extracts, Other biological fluids, Platform: Microplate reader 560 nm
Sample Types	Urine, Serum, Plasma, Tissue extracts, Other biological fluids
Expiration Date	Please enquire.



Representative
Standard
Curve