

## Product Datasheet

# Rabbit Human IgG F (ab')<sub>2</sub> antibody (FITC) (orb1473583)

<b>Description</b>	Rabbit polyclonal antibody against Rabbit Human IgG F (ab') <sub>2</sub> conjugated to
<b>Species/Host</b>	Rabbit
<b>Reactivity</b>	Human
<b>Conjugation</b>	FITC
<b>Tested Applications</b>	ELISA, FC, IF, IH
<b>Immunogen</b>	Human IgG F(ab') <sub>2</sub>
<b>Preservatives</b>	Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, containing 1% rAlbumin, 50% glycerol, 0.05% Sodium Azide
<b>Form/Appearance</b>	1 mg/ml. Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, containing 1% rAlbumin, 50% glycerol, 0.05% Sodium Azide
<b>Storage</b>	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
<b>Note</b>	For research use only
<b>Clonality</b>	Polyclonal
<b>Source</b>	Rabbit
<b>Dilution Range</b>	IH (1/100 - 1/1000), IF (1/100 - 1/1000), FC (1/1000 - 1/4000)
<b>Expiration Date</b>	12 months from date of receipt.