

## Product Datasheet

# Goat Human IgG F (ab')<sub>2</sub> antibody (FITC) (orb1473584)

|                            |  |
|----------------------------|--|
| <b>Description</b>         | Goat polyclonal antibody against Goat Human IgG F (ab') <sub>2</sub> conjugated to                                   |
| <b>Species/Host</b>        | Goat   |
| <b>Reactivity</b>          | Human  |
| <b>Conjugation</b>         | FITC   |
| <b>Tested Applications</b> | ELISA, FC, IF, IH  |
| <b>Immunogen</b>           | Human IgG F(ab') <sub>2</sub>  |
| <b>Preservatives</b>       | Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, containing 1% rAlbumin, 50% glycerol, 0.05% Sodium Azide          |
| <b>Form/Appearance</b>     | 1 mg/ml. Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, containing 1% rAlbumin, 50% glycerol, 0.05% Sodium Azide |
| <b>Storage</b>             | Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.                     |
| <b>Note</b>                | For research use only  |
| <b>Clonality</b>           | Polyclonal   |
| <b>Source</b>              | Goat   |
| <b>Dilution Range</b>      | IH (1/100 - 1/1000), IF (1/100 - 1/1000), FC (1/1000 - 1/4000)   |
| <b>Expiration Date</b>     | 12 months from date of receipt.  |