



Product Datasheet

FAM120A antibody (orb1473651)


Description	Mouse monoclonal antibody to FAM120A
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, IH, WB
Immunogen	Recombinant fusion protein of human FAM120A. The exact sequence is proprietary.
Target	FAM120A
Preservatives	Mouse IgG1 kappa. Liquid in PBS, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Mouse IgG1 kappa. Liquid in PBS, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Monoclonal
Source	Mouse
Uniprot ID	Q9NZB2
Entrez	23196
Dilution Range	WB (1/500 - 1/2000), IH (1/50 - 1/200), IF/IC (1/10 - 1/50)
Expiration Date	12 months from date of receipt.



Western blot analysis of FAM120A express...



Immunohistochemical analysis of FAM120A ...



Immunofluorescent analysis of FAM120A st...