



## Product Datasheet

### CRMP-5 antibody (orb1473693)

<b>Description</b>	Mouse monoclonal antibody to CRMP-5
<b>Species/Host</b>	Mouse
<b>Reactivity</b>	Human, Mouse, Rat
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	IF, WB
<b>Immunogen</b>	Recombinant fusion protein of human CRMP-5. The exact sequence is proprietary.
<b>Target</b>	DPYSL5
<b>Preservatives</b>	Mouse IgG1 kappa. Liquid in PBS, pH 7.3, 30% glycerol, and 0.01% sodium azide.
<b>Form/Appearance</b>	Mouse IgG1 kappa. Liquid in PBS, pH 7.3, 30% glycerol, and 0.01% sodium azide.
<b>Storage</b>	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
<b>Note</b>	For research use only
<b>Clonality</b>	Monoclonal
<b>Source</b>	Mouse
<b>Uniprot ID</b>	<a href="#">Q9JHU0</a> , <a href="#">Q9BPU6</a> , <a href="#">Q9EQF6</a>
<b>Entrez</b>	<a href="#">56896</a> , <a href="#">100910964</a> , <a href="#">65254</a> , <a href="#">65208</a>
<b>Dilution Range</b>	WB (1/1000 - 1/2000), IF/IC (1/10 - 1/50)
<b>Expiration Date</b>	12 months from date of receipt.



Western blot analysis of CRMP-5 expressi...



Immunofluorescent analysis of CRMP-5 sta...