



Product Datasheet

p27Kip1 (Ab0 Antibody (orb14737))

Description	p27Kip1 (Ab0 Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	P46527
Expiration Date	12 months from date of receipt.



Western blot analysis of p27Kip1 antibody...



Immunohistochemical analysis of formalin...