



Product Datasheet

BANP antibody (orb1473730)

| | |
|----------------------------|---|
| Description | Mouse monoclonal antibody to BANP |
| Species/Host | Mouse |
| Reactivity | Human |
| Conjugation | Unconjugated |
| Tested Applications | IH, WB |
| Immunogen | Recombinant fusion protein of human BANP. The exact sequence is proprietary. |
| Target | BANP |
| Preservatives | Mouse IgG1 kappa. Liquid in PBS, pH 7.3, 30% glycerol, and 0.01% sodium azide. |
| Form/Appearance | Mouse IgG1 kappa. Liquid in PBS, pH 7.3, 30% glycerol, and 0.01% sodium azide. |
| Storage | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles. |
| Note | For research use only |
| Clonality | Monoclonal |
| Source | Mouse |
| Uniprot ID | Q8N9N5 |
| Entrez | 54971 |
| Dilution Range | WB (1/500 - 1/1000), IH (1/50 - 1/200) |
| Expiration Date | 12 months from date of receipt. |



Western blot analysis of BANP expression...



Immunohistochemical analysis of BANP sta...