



Product Datasheet

RAPTOR antibody (orb1473732)

Description	Mouse monoclonal antibody to RAPTOR
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IH, WB
Immunogen	Recombinant fusion protein of human RAPTOR. The exact sequence is proprietary.
Target	RPTOR
Preservatives	Mouse IgG1 kappa. Liquid in PBS, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Mouse IgG1 kappa. Liquid in PBS, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Monoclonal
Source	Mouse
Uniprot ID	Q8N122 , Q8K4Q0
Entrez	57521
Dilution Range	WB (1/1000 - 1/2000), IH (1/50 - 1/200)
Expiration Date	12 months from date of receipt.



Western blot analysis of RAPTOR expressi...



Immunohistochemical analysis of RAPTOR S...