


Product Datasheet

ACOX1 antibody (orb1473774)


Description	Mouse monoclonal antibody to ACOX1
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, IH, WB
Immunogen	Recombinant fusion protein of human ACOX1. The exact sequence is proprietary.
Target	ACOX1
Preservatives	Mouse IgG1. Liquid in PBS, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Mouse IgG1. Liquid in PBS, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Monoclonal
Source	Mouse
Uniprot ID	Q15067
Entrez	51
Dilution Range	WB (1/500 - 1/1000), IH (1/50 - 1/200), IF/IC (1/10 - 1/50)
Expiration Date	12 months from date of receipt.



Western blot analysis of ACOX1 expressio...



Immunohistochemical analysis of ACOX1 st...



Immunofluorescent analysis of ACOX1 stai...