

Product Datasheet

AKAP13 antibody (orb1473796)

Description	Mouse monoclonal antibody to AKAP13
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	KLH-conjugated synthetic peptide encompassing a sequence within the N-term region of human AKAP13. The exact sequence is proprietary.
Target	AKAP13
Preservatives	Mouse IgG1 kappa. Liquid in PBS, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Mouse IgG1 kappa. Liquid in PBS, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Monoclonal
Source	Mouse
Uniprot ID	Q12802
Entrez	11214
Dilution Range	WB (1/500 - 1/4000)
Expiration Date	12 months from date of receipt.



Western blot analysis of AKAP13 expressi...