



Product Datasheet

RAB11A antibody (orb1473834)

Description	Mouse monoclonal antibody to RAB11A
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IH, WB
Immunogen	Recombinant fusion protein of human RAB11A. The exact sequence is proprietary.
Target	RAB11A
Preservatives	Mouse IgG1 kappa. Liquid in PBS, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Mouse IgG1 kappa. Liquid in PBS, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Monoclonal
Source	Mouse
Uniprot ID	P62492 , P62491 , P62494
Entrez	81830 , 8766 , 53869
Dilution Range	WB (1/1000 - 1/2000), IH (1/50 - 1/200)
Expiration Date	12 months from date of receipt.



Western blot analysis of RAB11A expressi...



Immunohistochemical analysis of RAB11A S...