



Product Datasheet

PPP2R1B antibody (orb1473920)

Description	Mouse monoclonal antibody to PPP2R1B
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, IF, IH, WB
Immunogen	Recombinant fusion protein of human PPP2R1B. The exact sequence is proprietary.
Target	PPP2R1B
Preservatives	Mouse IgG1 kappa. Liquid in PBS, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Mouse IgG1 kappa. Liquid in PBS, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Monoclonal
Source	Mouse
Uniprot ID	Q4QQT4 , Q7TNP2 , P30154
Entrez	5519 , 315648 , 73699
Dilution Range	WB (1/500 - 1/2000), IH (1/50 - 1/200), IF/IC (1/10 - 1/50), FC (1/10 - 1/50)
Expiration Date	12 months from date of receipt.




Western blot analysis of PPP2R1B express...



Immunohistochemical analysis of PPP2R1B

...



Immunofluorescent analysis of PPP2R1B st...