

Product Datasheet

HNF1 alpha antibody (orb1473957)

Description	Mouse monoclonal antibody to HNF1 alpha
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	KLH-conjugated synthetic peptide encompassing a sequence within the center region of human HNF1 alpha. The exact sequence is proprietary.
Target	HNF1A
Preservatives	Mouse IgG1. Supplied in crude ascites with 0.01% sodium azide.
Form/Appearance	Mouse IgG1. Supplied in crude ascites with 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Monoclonal
Source	Mouse
Uniprot ID	P20823
Entrez	6927
Dilution Range	WB (1/500 - 1/1000)
Expiration Date	12 months from date of receipt.



Western blot analysis of HNF1 alpha expr...