



Product Datasheet

CD107b antibody (orb1474002)

Description	Mouse monoclonal antibody to CD107b
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	Recombinant fusion protein of human CD107b. The exact sequence is proprietary.
Target	LAMP2
Preservatives	Mouse IgG1 kappa. Liquid in PBS, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Mouse IgG1 kappa. Liquid in PBS, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Monoclonal
Source	Mouse
Uniprot ID	P13473
Entrez	3920
Dilution Range	WB (1/1000 - 1/2000)
Expiration Date	12 months from date of receipt.



Western blot analysis of CD107b expressi...



Immunofluorescent analysis of CD107b sta...