

Product Datasheet

C-Reactive Protein Antibody (orb1474102)

Catalog Number	orb1474102
Category	Antibodies
Description	The C-Reactive Protein Antibody is suitable for WB. It is a Monoclonal, Unconjugated antibody which raised against KLH-conjugated synthetic peptide encompassing a sequence within the N-term region of human C-Reactive Protein. The exact sequence is proprietary. Purification: This antibody is purified through a protein G column.
Clonality	Monoclonal
Species/Host	Mouse
Conjugation	Unconjugated
Reactivity	Human
Form/Appearance	Mouse IgG3. Liquid in PBS, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Purification	This antibody is purified through a protein G column.
Immunogen	KLH-conjugated synthetic peptide encompassing a sequence within the N-term region of human C-Reactive Protein. The exact sequence is proprietary.
UniProt ID	P02741
Tested applications	WB
Dilution range	WB (1/500 - 1/1000)
Antibody Type	Primary Antibody
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: +44 (0)1223 859353 | Fax: +1 (415) 651-8558

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558

Note For research use only

Entrez **1401**



Western blot analysis of C-Reactive Protein expression in NCIH460 (A) whole cell lysates. (Predicted band size: 25 kD; Observed band size: 28 kD)

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)