

Product Datasheet

CYP27B1 antibody (orb1474161)

| | |
|----------------------------|---|
| Description | Mouse monoclonal antibody to CYP27B1 |
| Species/Host | Mouse |
| Reactivity | Human |
| Conjugation | Unconjugated |
| Tested Applications | WB |
| Immunogen | KLH-conjugated synthetic peptide encompassing a sequence within the C-term region of human CYP27B1. The exact sequence is proprietary. |
| Target | CYP27B1 |
| Preservatives | Mouse IgG1 kappa. Liquid in PBS, pH 7.3, 30% glycerol, and 0.01% sodium azide. |
| Form/Appearance | Mouse IgG1 kappa. Liquid in PBS, pH 7.3, 30% glycerol, and 0.01% sodium azide. |
| Storage | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles. |
| Note | For research use only |
| Clonality | Monoclonal |
| Source | Mouse |
| Uniprot ID | O15528 |
| Entrez | 1594 |
| Dilution Range | WB (1/500 - 1/1000) |
| Expiration Date | 12 months from date of receipt. |



Western blot analysis of CYP27B1 express...