


Product Datasheet

S100-A6 antibody (orb1474265)

| | |
|----------------------------|---|
| Description | Mouse monoclonal antibody to S100-A6 |
| Species/Host | Mouse |
| Reactivity | Human |
| Conjugation | Unconjugated |
| Tested Applications | WB |
| Immunogen | Purified recombinant human S100A6 protein expressed in E.coli. |
| Target | S100A6 |
| Preservatives | Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide, pH 7.3. |
| Form/Appearance | Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide, pH 7.3. |
| Storage | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles. |
| Note | For research use only |
| Clonality | Monoclonal |
| Source | Mouse |
| Uniprot ID | P06703 |
| Entrez | 6277 |
| Dilution Range | WB (1/500 - 1/1000) |
| Expiration Date | 12 months from date of receipt. |



Western blot analysis of S100-A6 express...