


Product Datasheet

MKK6 antibody (orb1474270)

Description	Mouse monoclonal antibody to MKK6
Species/Host	Mouse
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	IP, WB
Immunogen	Purified recombinant human MKK6 protein fragments expressed in E.coli.
Target	MAP2K6
Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide, pH 7.3.
Form/Appearance	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide, pH 7.3.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Monoclonal
Source	Mouse
Uniprot ID	P52564
Entrez	5608
Dilution Range	WB (1/500 - 1/1000), IP (1/10 - 1/50)
Expiration Date	12 months from date of receipt.



Western blot analysis of MKK6 expression...