


Product Datasheet

TAK1 antibody (orb1474296)

Description	Mouse monoclonal antibody to TAK1
Species/Host	Mouse
Reactivity	Human, Mouse, Primate, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Purified recombinant human TAK1 protein fragments expressed in E.coli.
Target	MAP3K7
Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide, pH 7.3.
Form/Appearance	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide, pH 7.3.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Monoclonal
Source	Mouse
Uniprot ID	Q62073 , O43318 , P0C8E4
Entrez	6885 , 26409 , 100910771 , 313121
Dilution Range	WB (1/500 - 1/1000)
Expiration Date	12 months from date of receipt.



Western blot analysis of TAK1 expression...