


Product Datasheet

P2Y12 antibody (orb1474435)

Description	Rabbit monoclonal antibody to P2Y12
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IP, WB
Immunogen	A synthetic peptide of human P2Y12
Target	P2RY12
Preservatives	Liquid in 50mM Tris-Glycine (pH 7.4), 0.15M NaCl, 50% Glycerol, 0.01% Sodium azide and 0.05% rAlbumin.
Form/Appearance	Liquid in 50mM Tris-Glycine (pH 7.4), 0.15M NaCl, 50% Glycerol, 0.01% Sodium azide and 0.05% rAlbumin.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Monoclonal
Source	Rabbit
Uniprot ID	Q9H244
Entrez	64805
Dilution Range	WB (1/500 - 1/1000), IP (1/10 - 1/50)
Expiration Date	12 months from date of receipt.



Western blot analysis of P2Y12 expression...