


---

## Product Datasheet

**JNK1/2/3 (phospho-T183/Y185) antibody  
(orb1474463)**

**Description**

Rabbit monoclonal antibody to JNK1/2/3 (phospho-


**Species/Host**

Rabbit

**Reactivity**

Human, Rat

**Conjugation**

Unconjugated

**Tested Applications**

IP, WB

**Immunogen**

KLH-conjugated synthetic phosphopeptide corresponding to residues surrounding T183/Y185 of human JNK1/2/3 protein. The exact sequence is proprietary.

**Target**

MAPK8

**Preservatives**

Liquid in 50mM Tris-Glycine (pH 7.4), 0.15M NaCl, 50% Glycerol, 0.01% Sodium azide and 0.05% rAlbumin.

**Form/Appearance**

Liquid in 50mM Tris-Glycine (pH 7.4), 0.15M NaCl, 50% Glycerol, 0.01% Sodium azide and 0.05% rAlbumin.

**Storage**

Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.

**Note**

For research use only

**Clonality**

Monoclonal

**Source**

Rabbit

**Uniprot ID**
[P49185](#), [P45983](#)
**Entrez**
[116554](#), [5599](#)
**Dilution Range**

WB (1/500 - 1/1000), IP (1/10 - 1/50)

**Expiration Date**

12 months from date of receipt.