



Product Datasheet

GSK3β (Phospho-Ser9) Antibody (orb14745)


Description	GSK3β (Phospho-Ser9) Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	P49841
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of paraffin...



Western blot analysis of extracts from 2...



Immunofluorescence staining of methanol-
...