



Product Datasheet

BAG6 antibody (orb1474561)

Description	Rabbit monoclonal antibody to BAG6
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	A synthetic peptide of human BAT3
Target	BAG6
Preservatives	Liquid in 50mM Tris-Glycine (pH 7.4), 0.15M NaCl, 50% Glycerol, 0.01% Sodium azide and 0.05% rAlbumin.
Form/Appearance	Liquid in 50mM Tris-Glycine (pH 7.4), 0.15M NaCl, 50% Glycerol, 0.01% Sodium azide and 0.05% rAlbumin.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Monoclonal
Source	Rabbit
Uniprot ID	P46379 , Q9Z1R2 , Q6MG49
Entrez	224727 , 94342 , 7917
Dilution Range	WB (1/500 - 1/1000), IF/IC (1/50 - 1/100)
Expiration Date	12 months from date of receipt.



Western blot analysis of BAG6 expression...



Immunofluorescent analysis of BAG6 stain...