



Product Datasheet

ABI2 Antibody (orb1474864)


Description	ABI2 Antibody
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, FC, ICC, IF, WB
Immunogen	E.coli-derived human ABI2 recombinant protein (Position: H151-D400).
Form/Appearance	Lyophilized
Concentration	Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.
Storage	At -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freezing and thawing.
Note	For research use only
Application notes	Tested Species: In-house tested species with positive results. Other applications have not been tested. Optimal dilutions should be determined by end users. Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	56 kDa
Uniprot ID	Q9NYB9
Expiration Date	12 months from date of receipt.



IF analysis of SNF2H/SMARCA5 and Tubulin...



Western blot analysis of ABI2 using anti-ABI2 antibody.



Flow Cytometry analysis of K562 cells using anti-ABI2 antibody.