

---

## Product Datasheet

**Desmoglein 2/DSG2 Antibody (monoclonal, 2B4D1)  
(orb1474876)**

**Description**

Desmoglein 2/DSG2 Antibody (monoclonal,


**Species/Host**

Mouse

**Reactivity**

Human

**Conjugation**

Unconjugated

**Tested Applications**

IHC, WB

 IHC  
analysis of  
Desmoglein  
2/DSG2  
using ...

**Immunogen**

E.coli-derived human Desmoglein 2/DSG2 recombinant protein (Position: L24-E1020).

**Form/Appearance**


Lyophilized

**Concentration**

Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.

**Storage**

At -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freezing and thawing.



 IHC  
analysis of  
Desmoglein  
2/DSG2  
using ...

**Note**

For research use only

**Application notes**

Tested Species: In-house tested species with positive results. Other applications have not been tested. Optimal dilutions should be determined by end users. Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.


 IHC  
analysis of  
Desmoglein  
2/DSG2  
using ...

**Isotype**

IgG1

**Clonality**

Monoclonal

**MW**

160 kDa

**Uniprot ID**
[Q14126](#)
**Expiration Date**

12 months from date of receipt.