



Product Datasheet

IRBIT/AHCYL1 Antibody (monoclonal, 2E7D9) (orb1474880)

Description	IRBIT/AHCYL1 Antibody (monoclonal, 2E7D9)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, WB
Immunogen	E.coli-derived human IRBIT/AHCYL1 recombinant protein (Position: E14-K57).
Form/Appearance	Lyophilized
Concentration	Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.
Storage	At -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freezing and thawing.
Note	For research use only
Application notes	Tested Species: In-house tested species with positive results. Other applications have not been tested. Optimal dilutions should be determined by end users. Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.
Isotype	IgG1
Clonality	Monoclonal
MW	59 kDa
Uniprot ID	O43865
Expiration Date	12 months from date of receipt.



Western blot analysis of IRBIT/AHCYL1 us...



Flow Cytometry analysis of K562 cells us...