

---

## Product Datasheet

**Ikaros/IKZF1 Antibody (monoclonal, 5F12H7)  
(orb1474882)**

**Description**

Ikaros/IKZF1 Antibody (monoclonal, 5F12H7)

**Species/Host**

Mouse

**Reactivity**

Human

**Conjugation**

Unconjugated

**Tested**

WB

**Applications**
**Immunogen**

A synthetic peptide corresponding to a sequence at the C-terminus of human Ikaros (428-459aa LKEEHRAYDLLRAASENSQDALRVVSTSGEQM), different from the related mouse sequence by five amino acids.

**Form/Appearance**

Lyophilized

**Concentration**

Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.

**Storage**

At -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freezing and thawing.

**Note**

For research use only

**Application notes**

Tested Species: In-house tested species with positive results. Other applications have not been tested. Optimal dilutions should be determined by end users. Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.

**Isotype**

IgG1

**Clonality**


Monoclonal

**MW**

55-65 kDa

**Uniprot ID**
[Q13422](#)
**Expiration Date**

12 months from date of receipt.



Western blot analysis of Ikaros/IKZF1 us...