

## Product Datasheet

# Rabbit Anti-Rat IgG (H+L) Secondary Antibody, Unconjugated (orb1474894)

<b>Description</b>	Rabbit Anti-Rat IgG (H+L) Secondary Antibody, Unconjugated
<b>Species/Host</b>	Rabbit
<b>Reactivity</b>	Rat
<b>Conjugation</b>	Unconjugated
<b>Immunogen</b>	Rat IgG (whole molecule).
<b>Storage</b>	Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.
<b>Note</b>	For research use only
<b>Clonality</b>	Polyclonal
<b>Expiration Date</b>	12 months from date of receipt.