

Product Datasheet

Rabbit Anti-Goat IgG (H+L) Secondary Antibody, Unconjugated (orb1474896)

Description	Rabbit Anti-Goat IgG (H+L) Secondary Antibody, Unconjugated
Species/Host	Rabbit
Reactivity	Goat
Conjugation	Unconjugated
Immunogen	Goat IgG (whole molecule).
Form/Appearance	lyophilized
Storage	Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Expiration Date	12 months from date of receipt.