

Product Datasheet

Goat Anti-Rabbit IgG (H+L) Secondary Antibody, Unconjugated (orb1474897)

Description	Goat Anti-Rabbit IgG (H+L) Secondary Antibody, Unconjugated
Species/Host	Goat
Reactivity	Rabbit
Conjugation	Unconjugated
Immunogen	Rabbit IgG (whole molecule).
Storage	Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
Expiration Date	12 months from date of receipt.