



Product Datasheet

Goat Anti-Rabbit IgG (H+L) Secondary Antibody, TRITC Conjugated (orb1474902)

Description	Goat Anti-Rabbit IgG (H+L) Secondary Antibody, TRITC Conjugated
Species/Host	Goat
Reactivity	Rabbit
Conjugation	TRITC
Tested Applications	FC, IF
Immunogen	Rabbit IgG (whole molecule).
Concentration	1mg/ml
Storage	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Note	For research use only
Clonality	Polyclonal
Expiration Date	12 months from date of receipt.



IF analysis of human mammary cancer tiss...



IF analysis of human intestinal cancer t...