


## Product Datasheet

# Goat Anti-Mouse IgG (H+L) Secondary Antibody, TRITC Conjugated (orb1474903)

<b>Description</b>	Goat Anti-Mouse IgG (H+L) Secondary Antibody, TRITC Conjugated
<b>Species/Host</b>	Goat
<b>Reactivity</b>	Mouse
<b>Conjugation</b>	TRITC
<b>Tested Applications</b>	FC, IF
<b>Immunogen</b>	Mouse IgG (whole molecule).
<b>Concentration</b>	1mg/ml
<b>Storage</b>	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
<b>Note</b>	For research use only
<b>Clonality</b>	Polyclonal
<b>Expiration Date</b>	12 months from date of receipt.



IHC analysis of rat heart using anti-Alp...