

Product Datasheet

HRP Conjugated AffiniPure Rabbit Anti-Rat IgG (H+L) (orb1474913)

Description	HRP Conjugated AffiniPure Rabbit Anti-Rat IgG (H+L)
Species/Host	Rabbit
Reactivity	Rat
Conjugation	HRP
Tested Applications	ELISA, WB
Immunogen	Rat IgG (whole molecule).
Form/Appearance	Concentrated, Liquid
Concentration	1mg/ml
Storage	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing.
Note	For research use only
Application notes	Diluted by neutral PBS or TBS. Use of sodium azide as a preservative will substantially inhibit the enzyme activity of horseradish peroxidase.
Isotype	IgG
Clonality	Polyclonal
Expiration Date	12 months from date of receipt.