

Product Datasheet

Rabbit Anti-Rat IgG (H+L) Secondary Antibody, FITC Conjugated (orb1474923)

Description	Rabbit Anti-Rat IgG (H+L) Secondary Antibody, FITC Conjugated
Species/Host	Rabbit
Reactivity	Rat
Conjugation	FITC
Tested Applications	FC, IF
Immunogen	Rat IgG (whole molecule).
Concentration	1mg/ml
Storage	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Note	For research use only
Clonality	Polyclonal
Expiration Date	12 months from date of receipt.