

Product Datasheet

Goat Anti-Rabbit IgG (H+L) Secondary Antibody, FITC Conjugated (orb1474924)


Description	FITC Conjugated Goat Anti-rabbit IgG (H+L) secondary antibody This FITC conjugated antibody is specific for rabbit IgG and shows no cross-reactivity with mouse/bovine IgG.
Species/Host	Goat
Reactivity	Rabbit
Conjugation	FITC
Tested Applications	FC, IF
Immunogen	Rabbit IgG (whole molecule).
Concentration	1mg/ml
Storage	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Note	For research use only
Clonality	Polyclonal
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)



IF analysis of rat brain tissues using rabbit anti-PKC Iota Antigen Affinity purified polyclonal antibody. FITC Conjugated Goat Anti-rabbit IgG (H+L) secondary antibody was used to detect the primary antibody.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)