

## Product Datasheet

# Mouse Anti-Rabbit IgG (H+L) Secondary Antibody, Biotin Conjugated (orb1474926)

<b>Description</b>	Mouse Anti-Rabbit IgG (H+L) Secondary Antibody, Biotin Conjugated
<b>Species/Host</b>	Mouse
<b>Reactivity</b>	Rabbit
<b>Conjugation</b>	Biotin
<b>Tested Applications</b>	ELISA, IHC
<b>Immunogen</b>	Rabbit IgG (whole molecule).
<b>Concentration</b>	1mg/ml
<b>Storage</b>	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing.
<b>Note</b>	For research use only
<b>Clonality</b>	Polyclonal
<b>Expiration Date</b>	12 months from date of receipt.