

Product Datasheet

Goat Anti-Mouse IgM Secondary Antibody, Biotin Conjugated (orb1474934)

Description	Goat Anti-Mouse IgM Secondary Antibody, Biotin Conjugated
Species/Host	Goat
Reactivity	Mouse
Conjugation	Biotin
Tested Applications	ELISA, IHC
Immunogen	Mouse IgM (whole molecule).
Concentration	1mg/ml
Storage	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing.
Note	For research use only
Clonality	Polyclonal
Expiration Date	12 months from date of receipt.