


Product Datasheet

Human NT-proBNP Protein (orb1474960)

Description	N-terminal pro-brain (or B-type) natriuretic peptide (NT-proBNP) is produced predominately by the c...
Reactivity	Human
Endotoxins	< 0.2 EU/ug
Tested Applications	ELISA, WB
Preservatives	Add deionized water to prepare a working stock solution of approximately 1 mg/mL and let the lyophilized pellet dissolve completely
Form/Appearance	Lyophilized at 1 mg/mL in PBS
Storage	Store lyophilized protein at -20°C. Aliquot reconstituted protein and store at -80°C. Avoid repeated freezing/thawing cycles
Tag	N-terminal 6xHis
Note	For research use only
Protein Sequence	MRGSHHHHHHGMASMTGGQMGRDLYDDDDKDRWGSHPLGSPGSASDLET SGLQEQRNH LQGKLSSELQVEQTSLEPLQESPRPTGVWKSREVATEGIRGHRKMVLYTLRAPR
Purity	> 95%
Source	E.coli
MW	12.6 KDa
Expiration Date	6 months from date of receipt.



BCA to determine the quantity of the pro...