

Product Datasheet

Human IL-17RA Protein (Biotin) (orb1475203)

Description	Biotinylated Recombinant Human Interleukin-17 receptor A is produced by our Mammalian expression sys...
Endotoxins	Less than 0.1 ng/μg (1 EU/μg) as determined by LAL test.
Form/Appearance	Lyophilized from a 0.2 um filtered solution of PBS,pH7.4.
Storage	Lyophilized protein should be stored at ≤ -20°C, stable for one year after receipt. Reconstituted protein solution can be stored at 2-8°C for 2-7 days. Aliquots of reconstituted samples are stable at ≤ -20°C for 3 months.
Tag	C-6His-Avi
Note	For research use only
Uniprot ID	Q96F46
Expiration Date	6 months from date of receipt.