

Product Datasheet

Human HER3 Protein (Biotin) (orb1475226)

Description	Biotinylated Recombinant Human Receptor Tyrosine-Protein Kinase ErbB-3 is produced by our Mammalian ...
Endotoxins	Less than 0.1 ng/μg (1 EU/μg) as determined by LAL test.
Form/Appearance	Lyophilized from a 0.2 um filtered solution of PBS, pH 7.4.
Storage	Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.
Tag	C-6His-Avi
Note	For research use only
Uniprot ID	P21860
Expiration Date	6 months from date of receipt.