

Product Datasheet

Human Siglec-10 Protein (Biotin) (orb1475231)

Description	Biotinylated Recombinant Human Sialic Acid-binding Ig-like Lectin 10 is produced by our Mammalian ex...
Endotoxins	Less than 0.1 ng/μg (1 EU/μg) as determined by LAL test.
Form/Appearance	Lyophilized from a 0.2 um filtered solution of PBS,pH7.4.
Storage	Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.
Tag	C-Fc-Avi
Note	For research use only
Uniprot ID	Q96LC7
Expiration Date	6 months from date of receipt.