

Product Datasheet

Human Siglec-2 Protein (Biotin) (orb1475242)

| | |
|------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Description | Biotinylated Recombinant Human B-cell Receptor CD22 is produced by our Mammalian expression system a... |
| Endotoxins | Less than 0.1 ng/μg (1 EU/μg) as determined by LAL test. |
| Form/Appearance | Lyophilized from a 0.2 um filtered solution of PBS,pH7.4. |
| Storage | Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months. |
| Tag | C-6His-Avi |
| Note | For research use only |
| Uniprot ID | P20273 |
| Expiration Date | 6 months from date of receipt. |